
HOFFMANN'S


PHOSPHAT CEMENT normalhärtend normal setting

CE 0482

MADE IN GERMANY

1

2

3

HOFFMANN'S CEMENT

normalhärtend

Zinkphosphat-Cement

de

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Zinc phosphate cement

en

Instruction for use for the dentist

Composition:

Powder contains: zinc oxide, magnesium oxide;

Liquid contains: o-phosphoric acid.

Areas of Application:

Cavity lining as basis for all filling materials (amalgam, composites etc.); cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal as well as metal ceramics (zircon, aluminium and lithium disilicate ceramics) and crowns and bridges on implants; cementation of orthodontic bands, retention pins and screws, cement build up, temporary fillings.

Dispensing Instructions:

Cementation consistency: 1.5 g powder : 1.0 g liquid

Cavity lining consistency: 2.1 g powder : 1.0 g liquid

Dispensing and testing climate for testing in accordance with DIN EN ISO 9917:

Mix quantities indicated at 23 ± 1°C and a relative humidity of 50 ± 10%.

The dispensing instructions also apply when mixing the liquid with HOFFMANN'S UNIVERSAL CEMENT.

Processing Recommendations:

A clean, dry glass plate at ca. 20°C (68 °F) and a non-discolouring spatula are recommended for mixing.

Mixing:

In order to achieve a homogenous consistency, divide the measured quantity of powder into 4 portions (½ ¼ ¼ ¼). Beginning with the smallest quantity, stroke the powder into the measured liquid within 90 seconds. The cavity lining consistency has been reached when on lifting the spatula, the peak drawn out forms a small hook and does not sink back into the paste (Fig. 1). The cementation consistency has been reached when the peak formed on lifting the spatula slowly falls back into the paste (Fig. 2 and 3). Discard remaining powder or impure liquid.

Mixing Time: 90 seconds

Processing Time from beginning of mixing (at 23 °C, 72°F):

For cementation: 3:00 minutes

For cavity linings: 2:30 minutes

Setting Time from end of mixing:

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For cavity linings: 3:30 – 5:30 minutes

Special Information:

- It is vital that opened pulp is covered.
- By deep cavity short-term irritation of the pulp could occur, however, this irritation will recede completely after a few minutes. For cavity lining we recommend to use Hoffmann's Copal Varnish (Ord.-Nr. 81901).
- To obtain an optimal pH-value, comply with the mixing time of 90 s.
- The cement should be applied to a preferably dry tooth surface.
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.
- Close containers immediately after removing material.

Package Sizes:

Powder: 100 g, Liquid: 40 ml

Set: 35 g powder and 15 ml liquid

HOFFMANN'S CEMENT normal setting

- Do not use after the expiration date
- Compliant with DIN EN ISO 9917
- Material safety data sheet (MSDS) available at www.hoffmann-dental.com

* 15 Colours/Shades:

- 01 (white), -02 (bluish white), -03 (yellowish white), -04 (light yellow), -05 (yellow), -06 (gold yellow), -07 (gold brown), -08 (pearl grey), -09 (grey), -10 (greenish grey), -11 (bluish grey), -12 (brown), -13 (greyish brown), -14 (light pink), -15 (pink)

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HOFFMANN DENTAL MANUFAKTUR GmbH