



**HOFFMANN'S**



# Welcome

to the new and fully revised edition of Hoffmann's complete catalogue in the eventful history of the company. The story began in 1892 with the invention of zinc phosphate cement by chemist and founder Dr. Otto Hoffmann.

Hoffmann's phosphate cement – like the ginkgo leaf on the cover of our catalogue – is a strong and durable classic. The gold standard among dental cements is an integral part of the range of materials used by practices all over the world with an eye on sustainability.

We deliberately avoid the use of nanofillers or the modification of cements with monomers and other plastics.

Our sustainable and reliable materials for conservative and restorative dentistry are manufactured in Berlin using craftsmanship, modern processes and natural raw materials. When developing new products, we increasingly rely on proven natural remedies such as apitherapy and ozone therapy.

Join us on a voyage of discovery. Be brave to try new ways. Because today's innovations are tomorrow's baseline.

Enjoy reading!









## THE MANUFACTORY

Hoffmann's products are made with precise craftsmanship and great attention to detail. In the 19<sup>th</sup> century the master and his stoker toiled for several days in order to ensure the successful firing of their cement. The furnace was fed with wood and coal and heated to a maximum temperature of 1400 °C. Ever since, the raw material mixture is filled into saggars per hand..Today, however, modern technologies are used and combustion takes place in a high-performance computer-controlled furnace. We are proud to say that our production processes are also CO<sub>2</sub>-neutral. These are the qualified people that are entrusted

with the processes of manufacturing and distribution. They are going the extra mile every day.

Our multicultural team of employees allows us to accept different perspectives and understand the different needs related to dental products. As an entrepreneurial family, we anticipate things in the long term. This includes flexible working hours for parents and grandparents as well as investment in education and training.

What unites us is the willingness to learn from nature and the belief in sustainability.





# Quality Spanning Three Generations

## Award-winning safety

Already in the 19<sup>th</sup> century the inventor of phosphate cement, Dr. Otto Hoffmann, formulated the requirements that are still valid today according to which high-quality dental products should be produced on an industrial scale in order to make them available to as many people as possible for professional dental provision.

His striving for quality was already awarded the silver medal at the FDI's fifth International Dental Congress in Berlin in 1909. The second generation consistently continued to develop these high quality standards. Dietbert Hoffmann was an active member in national and international standardization organizations for the generation and enhancement of material standards from 1967 onwards.

Today, all activities are pooled in the ISO/TC 106 Dentistry, in which Hoffmann Dental Manufaktur is involved. A quality management system, which represents the entire company, was introduced in 1995 by today's Managing Director Tobias Hoffmann. Close monitoring of the entire value creation chain is ensured through strict

control of raw materials and end products, as well as regular in-process controls. All testing is carried out in the company's in-house laboratory.

In order to make sure that if it says Hoffmann's on it, there is also Hoffmann's in it, Hoffmann Dental Manufaktur was the first branded company in the dental industry to introduce product packaging with tamper-proof features. The combination of a traditional look and feel and modern standards was honoured with an award in the scope of the 9th International Corporate Design Award in 2005.

The employees' wealth of experience and the quality awareness of every individual are the motor for the continuous further development of the high standards of quality of this family enterprise.

Hoffmann's products are CE compliant, FDA approved and, since 2014, also certified as being Halal and Kosher. However, for Hoffmann Dental Manufaktur the adherence to statutory regulations and material standards is merely one of the minimum requirements.



DR. OTTO HOFFMANN  
1854 - 1938



# Inventor of phosphate cement

Chemist, company founder, bearded creator

Today, the portrait of the ingenious inventor can be found in dental practices all over the world. But who was Otto Hoffmann?

Otto Hoffmann was born in a village at the foot of the Harz mountains in 1854 as the son of a farmer. When he was 10 years old, the family moved to Halle, where Otto attended the boarding school of the Franckesche Stiftungen after the death of his father. Despite financial hardships, he studied in Berlin and Leipzig and received his doctorate in 1878 in Würzburg. As an independent chemist, he settled in Berlin and, in addition to dyes and many other useful preparations, invented a ground-breaking product for dentistry: zinc phosphate cement. From 1892 onwards, the innovative cement spread from Berlin all over the world.

Otto Hoffmann also liked to travel. In 1900 he visited the World Exhibition in Paris, and Spits-

bergen in 1903. But he was not only attracted by far-away places. In fact, ever increasing heights also fascinated him.

As a member of the Berliner Verein für Luftschiffahrt (Berlin Airship Association), he was an active balloonist. And countless mountain tours took him to the highest peaks of the Alps.

Perhaps it was the beard that made him so irresistible? At a ball of the Berlin Alpine Club in 1904 Otto met the love of his life, Lucie, who was 28 years younger than him. He then founded a family with her and built a house in Grunewald.

However, his work remained Otto Hoffmann's elixir of life, and he continued his vocation as a chemist until he was 83 years old.

His work as a pioneer in the dental industry lives on today, not least in the products bearing his name: Hoffmann's.



Read the series  
**HOFFMANN'S Tales**  
the true story of the inventor.  
Free ordering see page 72.



# The History of Zinc Phosphate Cement

Industrial revolution in dentistry

*“Let me sitfast in the inside of the tooth and give me its flesh as a home. From the tooth I will suck its blood and from the gums I will suck the dental pulp. ”*



This was written on an ancient Babylonian clay panel. At the time (approx. 1880 BC), “tooth worms” were seen as the cause of tooth aches. It was however still to take several thousand years until one would find the real causes of toothaches and begin taking scientific measures against these.

The revolutionising events of the industrial revolution, which pushed on rapid developments in medicine and dentistry, had the consequence that more and more professional solutions were sought after in dental provision. Eventually a virtually frenetic search began for suitable filling materials, starting in the mid 50s of the 19th century. Many a composition was already viewed rather critically by contemporaries and - from today's point of view - can only be seen as a hazardous concoction of alchemist ideas. Contemporary sources reported how lime hydrate was extracted from the ground dental enamel of carnivorous animals and made into a “paste” with phosphoric acid. The

following was noted about a mass called “dental kit”, consisting of 13 parts of finely ground quicklime and 12 parts of phosphoric acid: “The mixture can only be used within a time frame of 1-2 minutes and cannot be stirred due to its disintegrative properties.”

The first useful material was the “Cement” recommended by the royal Dresden dentist Augustin Rostaing from 1858 onwards and manufactured by his son, the chemist Charles Sylvester Rostaing - also known as Rostaing Dentinogenes or Rostaing Cement Plombage. The formula for these dentinogenes was kept absolutely secret by its two inventors. One could receive 12 ounces of the substance after one had deposited 20 Pounds Sterling or 155 Taler - a horrendous amount in those days - with a Dresden bank.


That is why dentists of the time complained not only about the price of the compound but also about its poor availability. When father and son died shortly after another, they took the produc-

tion secrets with them to their grave. The chemical base substances were actually known (zinc oxide and phosphoric acid), whereas the production in a furnace was an impossible task for the individual dentist.

The manufacture of useful cements was however the necessary prerequisite for the development of inlays.

Based on the basic ideas conceived by Rostaring, Dr. Otto Hoffmann was the first person to develop a high-quality dental cement after countless trials, and to produce it reliably and economically on an industrial scale. He founded the Harvard Dental Manufacturing Company in 1892 together with the dentist Dr. Robert Richter, who was the

person who incited him into thinking about researching such a material in the first place. After leading dentists in America, Germany and Austria had deemed the cement to be good, the first 100 portions of Otto Hoffmann zinc phosphate cement were officially traded on 12th February 1892. Up to the First World War, phosphate cement took on a world wide monopoly position and remained an important material in dentistry for generations. Hardly any other material has an equally long or successful clinical record in dental practice. Since the end of certain trademark disputes in 2005, our very own zinc phosphate cement is carrying the name of its inventor: Hoffmann's.



KAISERLICHES



PATENTAMT.

## PATENTSCHRIFT

— № 175382 —

KLASSE 30 A. GRUPPE 12.

AUSGEBEN DEN 11. SEPTEMBER 1906.

DR. OTTO HOFFMANN IN BERLIN.

Verfahren zur Herstellung eines Zahnzements.

Patentiert im Deutschen Reich vom 1. Juni 1904 ab.

Die meisten der bisher gebräuchlichen, aus Zinkoxyd und Phosphorsäure hergestellten Zahnzemente haben den Nachteil, in den Mundsäuren nicht völlig unlöslich zu sein, so daß sie von diesen allmählich angegriffen und gelöst werden. Ein weiterer Nachteil dieser Zahnzemente ist ihr durch ihre völlige Undurchsichtigkeit bedingtes totes Aussehen, durch das sie sich von den natürlichen Zähnen wesentlich abheben.

Es ist nun gelungen, einen von diesen Mängeln durchaus freien Zahnzement herzustellen, der sich durch völlige Unlöslichkeit in organischen Säuren auszeichnet und infolge seiner durchscheinenden Beschaffenheit dem Zahnschmelz des natürlichen Zahnes bei richtiger Färbung täuschend ähnlich ist. Das Verfahren besteht darin, daß man basisches oder neutrales Aluminiumphosphat, Aluminiumpyrophosphat oder Aluminiumborat, Stoffe, welche sich bekanntlich in organischen Säuren nicht lösen, mit einer Lösung von Ortho-, Meta- oder Pyrophosphorsäure

Der nach vorstehendem Verfahren bereitete Zahnzement ist gegen die Mundsäuren und sonstigen Einflüsse äußerst widerstandsfähig, hat infolge seines hohen Gehalts an Tonerde eine durchscheinende Beschaffenheit und hebt sich daher bei geeigneter Färbung von dem natürlichen Zahnschmelz nicht ab.

Beispielsweise löst man in 81 g Phosphorsäure vom spezifischen Gewicht 1,7 13 g Tonerdehydrat unter Zusatz von Wasser, filtriert und dampft auf 110 g ein; andererseits erhitzt man neutrales Aluminiumorthophosphat zu lebhafter Rotglut, pulvert es nach dem Erkalten und verrührt dieses Pulver für sich oder in Mischung mit etwa 30 Prozent von bei etwa 100° getrocknetem Tonerdehydrat mit obiger Flüssigkeit zu einer plastischen Masse. Oder man entwässert basisches Aluminiumphosphat ( $2 Al_2 O_3 \cdot P_2 O_5$ ) durch Erhitzen auf Rotglut, pulvert es sodann und verknetet es mit der wie angegeben bereiteten Phosphorsäure enthaltenden Flüssigkeit zu einer plastischen Masse.





# Copper

A trace element that really cleans up



Copper was the first metal that humans learned to form, approximately 10,000 years ago in the Neolithic Age. The word itself is derived from the Latin cuprum (aes cuprium), and refers to the island of Cyprus as a centre of copper mining in the ancient world.

The germicidal (bactericidal) effect of copper was known in the ancient world. The ancient Egyptians used copper for sterilizing drinking water, and Hippocrates of Kos described the treatment of open wounds with copper in his writings.

Today, the bactericidal effect of copper can not only be proven but also explained. As with silver, it has to do with the so-called oligodynamic effect: Copper ions cause a disruption of the bacterial metabolism and a degeneration of the DNA. And in fact, copper, unlike silver, even has the property of regenerating itself, therefore providing a permanently bactericidal effect.

Furthermore, copper is also one of the essential trace elements required by the human organism. According to the German Society for Nutrition, the recommended daily amount is between 1 and 1.5 milligrams. For example, ten grams of chocolate would cover the daily needs.

## Bactericidal protective and healing cement

Copper-containing alloys have been used successfully in dentistry for decades. They provide long-lasting restorative care and very good protection against secondary caries.

Phosphate cements with the addition of copper have been protecting against secondary caries for more than a century. In 1977 Schmalz proved that copper cements have a sustained bactericidal effect, which significantly exceeds the effect of corresponding cements without copper addition.

The antibacterial effect is based on the release of copper ions. The continuous release of copper ions leads to a significant reduction of microbial activity in the vicinity of the treated tooth, and a lasting positive effect.

Milk teeth are the traditional areas of application for copper ion cements. By using copal varnish the protective effect can be doubled – mechanically and chemically.

A relatively new area of application is the cementation of crowns and bridges on implants to protect against peri-implantitis.



Read the series  
**HOFFMANN'S Tales**  
the exciting story in Africa.  
Free ordering see page 72.





Agathis Australis



# Copal

## Young amber to protect the pulp

The name “copal” has become something of a collective term for fossil and recent resins of varying characteristics and botanical origins. Copalli is a loan word from Nahuatl, the ancient Aztec language, with which the Native Americans described smoked resin, which was sacrificed to Gods on their altars. Optically, copal is quite similar to amber and can sometimes contain trapped insects or other small creatures, which makes it therefore a popular collector's item. Copal can be almost colourless and transparent or take on a lush yellow to red-brown tone.

Fossil copals can be found in almost all parts of the world, in deposits a few meters under the ground. They can be thousands of years old, but are still rather young in comparison to amber, which can be several million years old. These natural resins originate from the damaged bark of conifers or leguminous (pulse) plants. Copals are usually described according to their origin, such as Zanzibar Copal or Manila Copal, and only sometimes according to the plant from which they originate, such as Kauri Copal, named after the New Zealand kauri fruit (*agathis australis*). East African copals are said to be the ones with the highest quality, as their characteristics come closest to those of amber. Fossil resins, also

known as mature copals, are hard, hardly soluble and difficult to melt. Recently extracted copals, such as tree copal, which is directly harvested from trees as fresh resin, are much more important for industrial processing.

At the beginning of the 20<sup>th</sup> Century, thousands of tons of copal were transported to Europe for the production of varnishes, where the hardest copals were used to make the most resistant varnishes. Today, however, mainly synthetic materials are used in the varnish industry. Copals, like natural resins in general, are today only used in cases where old techniques are implemented (especially in painting) or where environmentally friendly materials are needed (such as linoleum manufacturing).

The copals used by Hoffmann Dental Manufaktur for the varnish with the same name, are recently obtained copals from Indonesia. Copal is liquefied with ethanol and receives a lighter or darker colour according to its colour composition, similar to a Single Malt Whisky.

Copal is used for the production of thermoplastic impression compounds, and gives the products their typical shine and firmness. Copal varnish particularly serves as a cavity varnish to protect the pulp.



Propolis



# Propolis

Antibiotic of bees



Propolis is used by bees to protect the bee colony from disease. They use it to seal openings and to isolate foreign matter that has been introduced into the hive. The name of the honey-yellow, sweet-smelling bee resin comes from the Greek words for «before» and «city», thus describing the protective function that propolis has for the beehive.

For the most part propolis consists of pollen balm and natural resin, which young worker bees gather from buds and sores of birch, beech, alder and poplar trees. In the afternoon, when the sun has softened the resin, the young bees take on the resin with their mandibles and transport it back to the hive sticking to their legs. The collected propolis can be harvested if the beekeeper places a foreign body in the hive. In order to protect their

habitat, the bees surround the foreign body with propolis. The foreign body and the propolis sticking to it can then be removed.

More than 300 different ingredients can be attributed to propolis, including iron, copper, selenium, essential oils, vitamin A and vitamin E. The antioxidating effect of propolis is primarily explained by the high number of esters and polyphenols, especially flavonoids. The latter are largely responsible for the immune-enhancing and anti-inflammatory properties of propolis and have been used in medicine for a long time.

Propolis is used in Hoffmann's Pulpine, Hoffmann's Pulpine NE and Hoffmann's Pulpine Mineral for caries profunda therapy and for direct pulp capping.

# Propolis – a natural remedy

“Propolis is therefore an indispensable part of my dental office.”



*Christine Albinger-Voigt*



## A miracle of nature

In summer, up to 60,000 bees live in a beehive in a very confined space. During the incubation period, the temperature is approx. 35 °C, together with very high humidity. These are basically ideal conditions for the spread of germs – but bees have found ways to protect themselves from diseases:

They coat all surfaces of their hive with propolis. In this way they keep their innermost living space, the nursery and the pantry, germ-free. They say beehives are the most sterile places in nature. Every day more than 100,000 „landings“ take place in a beehive. Bees bring not only valuable pollen and sweet nectar from the outside world, but also bacteria and fungal spores. Since the area around the flight hole is covered with propolis, the landing takes place on a decontamination mat, so to speak.

Unwanted intruders (even mice!) are killed by bees with bee venom. If they then cannot remove these from the hive because of their size, they are covered with a propolis layer. This prevents decay and rot, and the intruders are actually preserved for years.

While bacteria repeatedly develop resistance to antibiotics, there is still no resistance to the natural antibiotic of bees. The reason for this is that the composition of propolis varies. Depending on what is available in the immediate surroundings of the hive, the bees extract the propolis from the young buds and shoots of different trees and shrubs. According to which criteria bees collect propolis is still a secret that has not yet been deciphered!

In any case, their success proves them right: bees have managed to survive over 100 million years and have spread all over the planet. Similarly to the way propolis has been protecting the superorganism of bees for millions of years, it also helps the human immune system.

## Effects of propolis

Propolis is therefore a valuable gift from bees – with a wide range of effects and applications. In South Korea, propolis even has its own research institute, which, interestingly enough, was a spin-off from the Institute for Nuclear Research. In the search for active substances that mitigate the effects of radioactive radiation, Korean nuclear physicians discovered propolis.

## ORAL AND DENTAL APPLICATION

Due to its antibacterial and anti-inflammatory properties, propolis is an ideal active ingredient for oral and dental care.

### Inflammations in mouth and throat

The disinfecting and anti-inflammatory properties of propolis are also appreciated for the treatment and relief of sore throats. Tea with honey is a popular household remedy, although it is better not to dissolve the honey in hot tea, otherwise the good ingredients (especially heat-sensitive vitamins and enzymes) are damaged. A mixture of propolis solution and honey in a ratio of 1:1 alleviates painful symptoms and reduces inflammation. In addition, the affected mucous membranes are strengthened in their protective function against invading viruses and bacteria.

Another possibility is to slowly melt 5 to 10 drops of a propolis tincture (the liquid from Hoffmann's Pulpine NE is very well suited) on a sugar cube in the mouth. Tough cookies can even take Pulpine with Eugenol for a severe sore throat.

### Caries Profunda

Fortunately, Caries Profunda treatments are very rare in my practice. However, if they do occur, I have been using Pulpine NE for 15 years with

very good results. The preparation developed by Dr. Kaczmarek was recommended to me by my colleagues and I have been able to maintain the vitality of numerous deeply damaged teeth over the years. The positive benefit of propolis preparations in the case of Caries Profunda and open pulp is simply incomparable.

### Open pulp

Pulpine NE can also be applied as a thin layer to ground tooth stumps and thus prevents grinding traumas. The pulps remain vital and the teeth do not show a negative reaction. Literature often refers to allergic reactions. However, I have not been able to detect any allergic or whole-body reactions so far.

Every patient is informed by me about the positive effects of this apitherapy before treatment. I find such an explanation important simply because the positive expectations (placebo effect!) ideally also promote healing. For me, propolis is therefore an indispensable part of my dental office!

## BIOLOGICAL SANDWICH TECHNIQUE



### 1. Hoffmann's PULPINE NE

Propolis, calcium hydroxide + zinc  
for direct and indirect pulp capping

### 2. Hoffmann's PHOSPHATE CEMENT


With zinc as dentine replacement

### 3. Composite

you trust in



*Christine Albinger-Voigt is a dentist and alternative practitioner with her own dental practice in Bad Homburg. As a member of professional dental societies and chairperson of the International Society for Holistic Dentistry (GZM), she is particularly committed to the creation of interdisciplinary networks through lectures and qualification measures.*



Syzygium aromaticum



# Eugenol

## Herba Benedicta



*Hildegard von Bingen*

The pharmacologically active constituent of natural clove oil is called eugenol, which is extracted from cloves. Eugenol also occurs in other plants, such as allspice, cinnamon, laurel or the real avens root (*geum urbanum*) found in Northern Europe, however, not in similarly high concentrations.

The clove tree (*syzygium aromaticum*) is a type of plant belonging to the myrtle family (*myrtaceae*), which originally comes from South-East Asia. The oldest archaeological find of cloves in ceramic vessels in Syria dates back to 1700 BC. Caravans brought the essential oil to Europe via the Incense Road. In ancient times and the Middle Ages, clove oil was considered one of the treasures vaulted by monasteries and patricians. The analgesic, antiseptic properties of clove oil have been known for centuries and are, for exam-

ple, described in the writings of Galen of Pergamum (131 – 201 AD) and Hildegard von Bingen (1098 – 1179 AD). The analgesic effects against toothache, in particular, are also mentioned in the writings. Hildegard von Bingen even refers to the plant as a “herba benedicta”, a blessed herb.

Clove oil is therefore likely to be the oldest remedy for toothache, which is still today widespread in medicine chests and dental practices around the world. It is particularly often used in emergencies for temporary filling. It acts as a mild local anaesthetic, with analgesic, soothing and antibacterial properties.

Hoffmann's uses pure essential clove oil from Indonesia in its formulation for PROXI APEX, a root canal filling material as well as for Hoffmann's Pulpine, a biological capping material for the protection of pulp.



# Palaquium gutta

# Gutta-percha

## A sap-giving tree



Gutta-percha is closely related to natural rubber. The main provider of this natural polyterpene or polyisoprene is gutta-percha (*isonandra gutta*), which is found in the rain forests of Malaysia and Indonesia. Its name also comes from Malaysia: getah=rubber and percha=tree. The milky sap (latex) is usually extracted by cutting notches into the trees' bark. Procedures also exist in which gutta-percha is extracted from leaves and twigs by using specific solvents.

In the year 1843, the Englishman Montgomery came across this interesting raw material in Singapore and immediately recognized its beneficial properties. After he had brought back samples with him to England, the large-scale import to Europe began shortly afterwards. The impending extermination of the sap-giving trees was avoided by plantation cultivation. In the middle of the 19<sup>th</sup> century this natural product was favoured for the manufacture of golf balls, jewellery and decorative accessories. Due to its good insulating properties, gutta-percha was also used for decades as a coating for electric cables.

Gutta-percha already found its way into dentistry in 1847 and was used as a filling material

for the first time in the USA. Originally molten lead was used for this purpose. Gutta-percha replaced lead and was used experimentally in many areas of application. For today's daily work, two forms are still of relevance: gutta-percha sticks for root canal fillings, and sealing gutta-percha for temporary sealing of cavities.

Gutta-percha has exceptional thermal isolation properties, high biocompatibility and a clinically proven antibacterial effect (bacteria find it especially difficult to attach to its surface).

The stick material is heated before use until it becomes soft and kneadable at approx. 50 °C and can then be inserted into the cavity in a thermoplastic state. Gutta-percha can then be removed relatively quickly and without having to drill.

Hoffmann's gutta-percha sticks consist of natural gutta-percha, to which zinc oxide and bees wax are added. The gutta-percha used at Hoffmann's comes from Java, is almost white and distinguishes itself through its especially high degree of purity. In order to maintain the elasticity of the material, Hoffmann's gutta-percha is sold only in re-sealable tubes.



*Fucus vesiculosus* and *Alaria esculenta*



# Algae

## Power from the sea

Algae are the oldest plants on our planet. Researchers assume a global total of more than 100,000 different species of algae. The algae filter valuable trace elements and minerals out of the sea through osmosis. Their nutrient density (amino acids, mineral salts, trace elements and vitamins) is therefore unique in the plant world.

The health-promoting active ingredients are fucoidan, which is obtained mainly from *laminaria japonica*, and alginic acid. Fucoidan plays an important role in the regeneration of tissue cells

and has a strong anti-inflammatory potential. Alginic acid ensures that the cell walls of the algae are stabilized. The gel-like texture of alginic acid gives the algae intercellular strength and flexibility at the same time.

The sodium alginate used is obtained from brown algae, which are particularly rich in trace elements and fucoidan. The algae have optimal adhesive properties, ideal for use in MAREFIX denture adhesive powders.

### HOFFMANN'S MARE FIX

Hoffmann's MAREFIX is a taste-neutral denture adhesive made from pure natural alginic acid.



#### AREAS OF APPLICATION

- Denture stomatitis

#### ADVANTAGES

- very good compatibility
- tasteless
- stable, elastic adhesive medium
- protects the oral mucosa from pressure points and injuries

#### COMPOSITION

Alginate in high pharmaceutical quality

Hoffmann's MAREFIX	
Order No.	Package size
84068	10 g powder



Kerria lacca



Video on the subject:  
[www.hoffmann-dental.com](http://www.hoffmann-dental.com)



# Shellac


## The work of an illustrious louse

Shellac is a symbiotic product of a tree and an insect which displays natural thermoplastic properties. Shellac has been used since the 15<sup>th</sup> century (dissolved in alcohol) for the finishing of furniture and instruments. It became world-famous at the end of the 19<sup>th</sup> century when it was used to make shellac records. However, the triumph of shellac records is not only due to Emil Berliner, the inventor of the gramophone, but equally due to an insect of barely 0.5 mm length called *coccus iacca kerr*, the shellac louse.


The term lac comes from the Sanskrit word „lakh“ which stands for „infinitely many“. Countless shellac lice settle in huge colonies on young tree shoots, especially those of *ficus religiosa* and *ficus indica*. The females of the shellac louse create lac by processing the resin of the fresh twigs in their organism. They themselves are bit by bit covered by the resin completely and die, while 20-30 larvae will develop inside of them, which

then escape as insects through holes. What remains is called seed lac, which encrusts the twigs tightly and can be harvested twice a year. The crust is removed from the twigs or branches mechanically and hot filtered through cloth. In order to produce one kilogram of shellac, one requires the metabolic products of approx. 300,000 shellac lice. Depending on the variety, this raw material is available in various colour variations from dark red and brown to yellow.


Hoffmann Dental Manufaktur uses only high-grade shellac from India. Hoffmann's base plates are produced in two different forms for upper and lower jaws. Dental technicians use them as carrier plates for model work for the setting of acrylic teeth or as a basis material for bite registration and individual/custom impression trays. Contrary to the light-curing and autopolymerizing synthetic material often used in dentistry, shellac is physiologically harmless and biologically degradable.



*CHRISTIAN FRIEDRICH SCHÖNBEIN*  
1799 – 1868





*WERNER VON SIEMENS*  
1816 – 1892



*NIKOLA TESLA*  
1856 – 1943



*JOACHIM HÄNSLER*  
1908 – 1981



# Ozone

## Nature's most effective disinfectant

Ozone - a colorless gas - is a molecule consisting of three oxygen atoms. It is the high-energy, active form of oxygen and is therefore also called "active oxygen". By supplying energy in the form of UV or electrical radiation, the oxygen is converted into ozone.

In the atmosphere, ozone is created by electrical discharge during lightning or water discharges. The ozone is responsible for the fresh and clean smell of the air that can be inhaled after a thunderstorm or near waterfalls. Ozone cleans and purifies air pollutants only for a short time, however, because it quickly converts back to normal oxygen. At the same time, this is an important advantage, as it literally accomplishes its task in a flash.

Science owes its discovery to an extraordinary German physicist and chemist, Christian Friedrich Schönbein, who discovered the gas in 1839 and named it „ozone“ after the Greek word for „smelling“, „ózon“.

Due to its instability, ozone cannot be stored for long periods of time or purchased in cylinders like other industrial gases. Before it can be used, it must be produced on site. There are different methods of generation. The first industrial ozone generators were developed by Werner von Siemens and Nikola Tesla.

### OZONE IN MEDICINE

Ozone is one of the most effective disinfectants and works even in very low concentrations in a relatively short time. Because of its high oxidation potential, ozone oxidizes the cellular components of bacteria cell walls and penetrates the cell.

Once it has penetrated, it oxidizes all essential components, thus damaging the membrane and killing off the cell. In addition to this bactericidal and viricidal effect, ozone also has sporicidal, fungicidal and protozoicidal properties.

Thanks to these properties, ozone has a firm place in preventive medicine, especially when it comes to the quality assurance of drinking water. Soon after its discovery, the gas established itself in medicine and is now used for the treatment of wounds, ulcers, burns as well as for the promotion of tissue healing and regeneration.

Already in 1920, A. Wolff used ozone water in the field of dentistry. Based on that, E. A. Fisch used ozone water in 1934 for the treatment of inflammatory changes. Joachim Hänsler then developed devices for ozone oxygen production for the medical sector.

### OZONIZED OILS

Ozonized oils can even be applied by patients at home if the instructions of the treating dentist are strictly followed. The amount of ozone used is taken into account both in the dental practice and at home and adapted to the clinical situation and development. The long shelf life (when cooled) of ozonized vegetable oils as well as their effectiveness make them particularly suitable for use in the dental practice as well as at home.



*Ricinus communis* and *olea europaea*



“Let food be thy medicine  
and medicine be thy food.”

Hippocrates

# Olive Oil

## The oldest remedy in the world



Olive oil is obtained from the fruit of *Olea Europaea* and consists mostly of monounsaturated oleic acid, and to a lesser extent of linoleic acid and palmitic acid.

Mainly two polyphenol antioxidants are responsible for the health-promoting benefits: oleuropein and oleocanthal. Thanks to oleuropein the olive tree is able to make free radicals harmless and counteract the pathogenesis of diseases. Oleuropein stimulates the immune system and slows the aging process of the olive tree. Oleocanthal is an analgesic aldehyde. This blocks the produc-

tion of neurotransmitters, which are responsible for inflammatory reactions.

The polyphenols in the olive oil, like other antioxidants, have an anti-inflammatory and health-promoting effect.

Additionally, vitamin E and squalane stimulate the cell renewal and prevent irritated skin with an impaired barrier function. Olive oil has particularly a positive effect on the gums and oral health. It helps to prevent bleeding of the gums and eliminates bacteria nested in the gingival sulcus.



# Castor Oil

## Healthy miracle



Ricinus oil is obtained from the seeds of *ricinus communis*, also known as the miracle tree. The medicinal plant, which is at home throughout the Mediterranean and wide-spread as an ornamental shrub, is already mentioned in the Old Testament: "God let the ricinus plant grow over Jonah for shady protection."

The seeds contain an extremely toxic ricin which is completely removed due to the pressing. Castor oil consists of 87% glycerides and ricinoleic acid. In addition, the oil has low levels of linoleic, palmitic and stearic acid. Castor

oil is especially known for its high content of tocopherol (vitamin E). In contrast to the olive oil ( $\alpha$ -tocopherol), however, castor oil is dominated by  $\gamma$ - and  $\delta$ -tocopherols. Both act as free radical scavengers and serve significantly as an antioxidant. The use of castor oil to combat pathogenic anaerobic bacteria in the oral flora has been proven. Hoffmann's PeriO<sub>3</sub> Oil provides ozonated castor and olive oil for the treatment of microbial diseases of the gums. The ozone-enriched oils have an antibacterial and anti-inflammatory effect. The restoration of a healthy oral flora is promoted.



The Gentle Alternative  
to Chlorhexidine



HOFFMANN'S PERIO<sub>3</sub> Oil



Natural gum thearapeutic based on ozonized olive oil and castor oil

AREAS OF APLICATION

- acute and chronic gingivitis
- periimplantitis
- stomatitis
- periodentitis
- cheilitis
- after professional tooth cleaning

ADVANTAGES

- natural preparation
- keeping the oral flora healthy
- plaque retardant
- 48 hours of oxygen release
- without antibiotics
- without added cortisone

HOW IT WORKS

Highly reactive ozonides are formed by the incorporation of ozone in natural oils. Unlike ozone gas, they are released slowly over a period of up to 48 hours. Anaerobic bacteria are selectively eliminated, and the oral flora regains its natural equilibrium. The use of chlorhexidine, but also of antibiotics, can be avoided or significantly reduced.



COMPOSITION

Ozone-enriched olive oil and castor oil

Hoffmann's Perio <sub>3</sub> Oil	
Order No.	Package size
84047	3 ml gel

# Bacteria-tight seal with PeriO<sub>3</sub> Oil



*"Some solutions are simple but ingenious"*


## PROF. DR. ÇETIN SEVÜK

*Prof. Dr. Çetin Sevük is an expert in prosthetics and implantology with decades of experience. In addition to his clinical work in his private clinic in Istanbul, he is active in several associations and academies and currently holds the position of President of Turkish Society for Prosthetics and Implantology.*

## Professor Sevük's advice for an osseointegration without complication

After placement of the implant screw, one drop of Hoffmann's PeriO<sub>3</sub> Oil is introduced into the screw gap (Figs. 1 and 2). A quantity of PeriO<sub>3</sub> Oil is then applied around the implant to seal the mucosa. Any sutured area is also covered with a thin layer of PeriO<sub>3</sub> Oil (Fig. 3).

After completion of the osseointegration period, Hoffmann's PeriO<sub>3</sub> Oil is applied to the implant screw gap and the surrounding tissue during insertion of the gingiva former and during the impression taking stages (Figs. 4 and 5). Some solutions are simple but ingenious.







**HOFFMANN'S GINGIVA FIT**



**AREAS OF APPLICATION**

- Support the treatment of inflammation and irritation of the oral mucosa and gums
- Support the treatment of periodontitis
- After the professional cleaning and after a curettage
- Inhibition of plaque formation by physical wetting of teeth and gums

**MECHANISM – HOW IT WORKS**

Gingiva Fit is produced by sparging ozone-oxygen through natural oils for several days thus forming ozonides. In the mouth the ozonides from the oil are slowly released over a period of 48 hours. Anaerobic bacteria are targeted and selectively reduced. At the same time the natural oils groom and stimulate the oral flora so that she can regain her natural healthy equilibrium.

Hoffmann's GINGIVA FIT	
Order No.	Package size
84048	1 ml gel

**ADVANTAGES**

- Natural product
- Conservation of oral flora
- Plaque inhibition
- Long duration of action
- Without antibiotics
- Without cortisone



**COMPOSITION**

Ozone-enriched olive oil and castor oil

# Biological periodontitis therapy at home



*"Our patients prefer Gingiva Fit to Chlorhexidine for home use"*

## PROF. DR. GLORIA SÁNCHEZ SABORIDO


Master of Implantology, Prosthetics and Periimplantology  
UNIVERSITY ALFONSO X EL SABIO  
Madrid, Spain

In all cases treated, a decrease in depth of penetration, inflammation, plaque index and bleeding index was observed. In the re-evaluations that we have carried out after scaling and root planing, we have collected impressions of patients regarding the use of Hoffmann's Gingiva Fit, all having previous experience with Chlorhexidine as a treatment adjuvant. In comparison with this antiseptic product they reported us following perceptions:

- absence of stinging / burning after gel application,
- absence of metallic taste or alteration of taste,
- absence of stains on the teeth,
- shorter treatment time than with Chlorhexidine.

Our patients prefer Gingiva Fit to Chlorhexidine for home use.

With practical resealable application cannula for home use.



Patients can apply Gingiva Fit to an interdental brush to access the interdental spaces.



## RENAISSANCE OF THE ZINC PHOSPHATE CEMENT



# 0%

BPA, monomeres,  
phenols, nanofillers



### HOFFMANN'S ZINC PHOSPHATE CEMENT – QUICK



#### AREAS OF APPLICATION

- Lining for all filling materials (amalgam, composites)
- Cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal as well as metal ceramic and all-ceramic (zirconium oxide, aluminium oxide and lithium disilicate ceramic)
- Cementation of implant-supported crowns and bridges
- Cementation of orthodontic bands
- Cementation of retention pins and screws
- Core build-ups
- Long-term temporary fillings

#### ADVANTAGES

- Dimensionally stable – no shrinkage and therefore no formation of marginal gaps
- Excellent micromechanical adhesion
- Wide range of applications from A like amalgam to Z like zirconia
- Very easy removal of excess material
- Extremely low pulp toxicity
- Hypoallergenic
- More than 127 years of clinical experience

#### COMPOSITION

Powder: zinc oxide, magnesium oxide

Cement liquid: o-phosphoric acid

#### PHOSPHATE CEMENT - quick

Order No.	Package size	Colour
80012-(Colour)	1 x 100g powder	01-15
8001300	1 x 40 ml liquid	---
80015-(Colour)	1 x 35 g powder 1 x 15 ml liquid	01-15

#### PHOSPHATE CEMENT - normal

Order No.	Package size	Colour
80022-(Colour)	1 x 100g powder	01-15
8002300	1 x 40 ml liquid	---
80025-(Colour)	1 x 35 g powder 1 x 15 ml liquid	01-15



HOFFMANN'S ZINC PHOSPHATE CEMENT – NORMAL



OVERVIEW OF TIMES

Times (min)

Areas of application	Phosphate Cement quick		
	Mixing time	Working time	Setting time
Luting	1:30	2:30	3:00 – 5:30
Lining	1:30	2:00	2:30 – 4:00

Areas of application	Phosphate Cement normal		
	Mixing time	Working time	Setting time
Luting	1:30	3:00	5:00 – 7:30
Lining	1:30	2:30	3:30 – 5:30

white	01
bluish white	02
yellowish	03
white	04
light yellow	04
yellow	05
gold yellow	06
gold brown	07
pearl grey	08
grey	09
greenish grey	10
bluish grey	11
brown	12
greyish brown	13
light pink	14
pink	15

Zn

ZINC

Hoffmann's zinc-phosphate-based cements replace lost dentin and additionally protect the pulp from thermal and bio-chemical stress.

- Remineralizing effect
- Demonstrably better than fluoride
- Protection against free radicals





# Mixing technique

A crown or bridge luted with Hoffmann's phosphate cement will be able to stay in a patient's mouth for decades. However, ideal success is only achieved if the cement is processed carefully. The two components, consisting of powder and liquid, are mixed together per hand.

This gives you a number of advantages, as the working time, amount and consistency can be individually adapted to the task at hand.

Hoffmann's provides two different dosage forms: the classic version in bottles for free-hand dispensing, and pre-dosed as READY2MIX in stick packs and fluid tubes.

Due to the exothermic reaction during the mixing procedure, the powder is always applied to the liquid in portions.

## Classic mixing by hand

For classic mixing, one divides the powder into four different portions ( $1/8$ ,  $1/8$ ,  $1/4$ ,  $1/2$ ). Starting with the smallest amount, the portions are mixed into the liquid one after the other with the spatula. The mass should be smoothed out across the slab every now and again, in order to allow the reaction warmth to escape well. A total

of 90 seconds is needed for the mixing process.

The most difficult but also the most important part is to determine the amount of powder which is to be put into the liquid. As a rule of thumb, a surplus of liquid should be avoided and the mixture should contain ample powder.

## Consistency control - Checkmark test method

The powder saturation is determined optically by means of the so-called checkmark test. To do so, a tip of material is pulled out of the mass with a spatula. If no tip can be pulled out or if the cement should even drip from the spatula, then powder will definitely have to be added.



### Lining (Pic. 1)

The correct lining consistency is achieved if the formed tip can be bent back into a hook, like a checkmark, and does not sink back into the mass.



### Luting

The correct luting consistency is achieved if the pulled out tip slowly sinks back into the mass. (Pic. 3) If it sinks back too fast, the mixture can

be adapted by adding a little more powder. If it does not sink back, but remains on the mass, (Pic. 2) then the consistency is too thick. In such a case, one should start mixing anew, as no more liquid should be added after mixing.

### Note

The cement should not be mixed any thinner for the luting of complicated cases that might require more time, as this might cause it to lose its stability. In such cases, instead of quick hardening cement, one can use normal hardening cement which provides a longer working time.



Pic. 1



Pic. 2



Pic. 3

# Non-toxic bonding


Powder and liquid of the phosphate cement essentially consist of zinc oxide and phosphoric acid. Zinc oxide is a component of wound creams for babies, phosphoric acid can be found in the most popular soft drink in the world: Coke.

The chemical components of phosphate cement are purely inorganic. The individual components are very well tolerated, which is also true for the end product. Zinc oxide and phosphoric acid react with each other and form a crystal, which is called hopeite:



The more powder is kneaded in, the more hopeite crystals are formed and the more favourable the cement's properties will become. The actual amount of hopeite crystals also has an effect on the solubility of the cement.

*The ancient Romans knew that only sand from rivers was suitable for the construction of houses. The least suitable is desert sand, the grains of which have been rounded through the strong winds.*



REM x 1000 magnification (Hoffmann's Cement, powder)

## MICROMECHANICAL BONDING

The adhesion of the phosphate cement occurs micromechanically. On the one hand, it is important that the tooth is prepared accordingly. Literature generally states an optimum preparation angle of 6°. On the other hand, the restoration should have a certain surface roughness.

An optimal particle size distribution ensures the best micromechanical properties of the cement. The distribution ideally follows a normal distribution with smaller and larger particles in the range of up to 25 microns. Additionally the surface roughness of the cement particles is another factor. Grinded particles without any roughness will not provide adhesion.

## ANGULAR POWDER PARTICLES

Hoffmann Dental Manufaktur already developed special grinder technology decades ago. The cements are ground to the right particle size in a three-stage grinding process, and under a microscope one can see that there are angular as well as sharp-edged particles.

## ACIDIC PAIN?

With zinc phosphate cements, so-called acid pain can occur with sensitive patients. This is caused by the acid effect in the direct vicinity of the pulp through nerve irritation. It has nothing to do with the warmth caused by the cement's reaction and immediately dissipates with the curing process.

## CAN BE AVOIDED

Acid pain can already be avoided by not reducing the recommended amount of powder and mixing it according to directions. In addition, the cavity can be prepared with cavity protection varnish, for example Hoffmann's Copal Varnish.





0%

BPA, monomeres,  
phenols, nanofillers



01	white
07	golden braun
10	greenish grey
11	blueish grey
15	pink



**HOFFMANN'S LEONARDO HARMONIC SHADES**



**UNIVERSAL GENIUS**

-  Aesthetic bonding of zirconia and lithium-disilicate
-  Dentin replacement / lining
- Semi-permanent aesthetic filling

TECHINCAL DATA		
	Hoffmann ´s	ISO standart
Setting time:	5:00 - 7:30 min	2:30 - 8:00 min.
Compressive strength:	90 - 130 MPa	50 Mpa
Film thickness:	< 20 µm	< 25 µm

HARMONIC SHADES / SET		
Order No.	Package size	Colour
82200	1 x 100 g, powder	01, 07,
	4 x 30 g, powder	10, 11,15
	1 x 40 ml, cement liquid	---
	1 x 40 ml, test fluid	---

**ADVANTAGES**

- Low allergy potential due to
  - purely inorganic chemistry
  - mineral food colours
- Optimum granulometry ensures
  - low film thickness
  - permanently stable micromechanical adhesion
- Fitting with test fluid
  - 100% colourfast
  - food quality
- Economical and ecological
  - no opaquer
  - no correction firing necessary

**COMPOSITION**

Powder: zinc oxide, magnesium oxide  
Cement liquid: o-phosphoric acid  
Test fluid: propane-1, 2, 3-Triol

## INDIVIDUAL SHADE CEMENTATION

Zinc phosphate cement has a decisive advantage especially for differently coloured tooth stumps, metallic post abutments and implant abutments – due to its opacity, it ensures an even priming and makes the translucent ceramics appear more even.

The set consists of one basic shade 01 and four strong colour shades in the direction of yellow, green, grey and red for individual shade adjustment. The cement shades are mixed individually. The colour effect of the various cements can also be determined using a digital tooth shade determination device.

## IMITATION OF VITAL PULP

“Pink trick” is the expression given to cementing with a pink cement. With this trick, pulp-dead teeth, whose grey tooth stumps cloud the overall impression of the restoration, are given a lively appearance, similar to pulp with blood circulation.



Basic colour 01  
+ shade 15



Basic colour 01

HARMONIC SHADES / REFILLS		
Order No.	Package size	Colour
82201	1 x 100 g powder	01 white
82507	1 x 30 g, powder	07 golden braun
82510	1 x 30 g, powder	10 greenish grey
82511	1 x 30 g, powder	11 blueish grey
82515	1 x 30 g, powder	15 pink
82300	1 x 40 ml, cement liquid	---
88800	1 x 40 ml, test fluid	---

Hoffmann's shades make it possible to modify the colour effect of the restoration „chair-side“. The test fluid included in the set can be used to test the optimum colour effect before the crowns are definitively cemented. This is particularly important in sensitive areas, e.g. in anterior tooth restoration. The cement mass with the test fluid does not set and can simply be rinsed off with water.

## MATCHING OF DIFFERENTLY SHADED TOOTH STUMPS

Lightening (base shade 01 / lightened shades with 01)



Darkening (shades 07, 10, 11)



Shade correction of the restoration  
(all shades and mixed shades)





# Hypoallergenic cementation for metal-free ceramics

Hoffmann's Harmonic shades  
Hoffmann's PHOSPHAT CEMENT normal setting  
Hoffmann's PHOSPHAT CEMENT quick setting  
Hoffmann's READY2MIX  
Hoffmann's READY2PROTECT Copperioncement



From an environmental medical point of view, they are the best choice for luting crowns and bridges made of zirconium oxide, aluminum oxide and lithium disilicat ceramic. The following table lists some well-known manufacturers and systems that are suitable for classical cementation with micro-mechanical adhesion. We update the list regularly, but make no claim to completeness.

Company	System	Name
3M ESPE AG	Lava <sup>TM</sup>	Lava <sup>TM</sup> Zirkonoxid Lava <sup>TM</sup> Plus hochtransluzentes Zirkonoxid
ACF GmbH	ZirLuna®	ZirLuna® ZirLuna® Sensitive
Amann Girrbach AG	ceramill®	ceramill® zolid fx classic ceramill® zolid fx preshades ceramill® zolid fx multilayer ceramill® zolid classic ceramill® zolid preshades ceramill® zi
Bien-Air GmbH	DC	DC-Zirkon® DC-Zirkon®Col DC-Leolux® DC-Procure® DC Shrink®
Bionah GmbH	BionZ	BionZ Crystal BionZ Diamond BionZ Glass BionZ Coral
CAD Esthetics AB	-	Denzir® Zirconia Magnum Zirconia Plus
DeguDent GmbH	Cercon®	Cercon®ht
Dental Direkt	DD Bio	DD Bio ZW iso (Color) DD Bio ZX², ZS, ZK DD Bio Cube X2 DD Bio ZA
Diadem SAS	-	Diazir

DOCERAM Medical Ceramics GmbH	Nacera®	Nacera® Pearl Multi-Shade Nacera® Pearl Shaded 16+2 Nacera® Standard 98 Nacera® Zirkonzahn* 95
Glidewell Laboratories	-	BruxZir® Solid Zirconia
Goldquadrat GmbH		Katana Zirconia ML Noritake Zirprime Quattro Disc Med
Ivoclar Vivadent AG	IPS e.max	IPS e.max Press IPS e.max ZirCAD IPS e.max CAD
Kavo Dental GmbH	-	Zirconia Soft
Kulzer	-	Cara Zirkon ZrST, ZrML, Zrtr Cara LiSi <sub>2</sub>
Kuraray Noritake Dental Inc.	KATANA™	KATANA™ Zirconia UTML / STML KATANA™ Zirconia ML KATANA™ Zirconia HT
Metoxit AG	-	Z-CAD®HD Z-CAD®HD HTL
Nobel Biocare AG	Procera®	Procera® Zirconia Procera® Alumina
Sagemax Bioceramics Inc.	NexxZr T	NexxZr T - Transluzentes Zirkon NexxZr S - Starkes Zirkon
Sirona Dental Systems GmbH	Cerec inCoris	Cerec Blocs C Cerec Blocs C PC inCoris ZI inCoris TZI inCoris TZI C inCoris ZI meso
Vita Zahnfabrik GmbH	In-Ceram	VITA SUPRINITY® VITA YZ HT VITA YZ T VITA In-Ceram YZ VITA In-Ceram AL VITA In-Ceram SPINELL VITA In-Ceram ZIRCONIA VITA In-Ceram ALUMINA
Wieland Dental+Technik GmbH & Co. KG	Zenostar®	Zenostar® MT Zenostar® T Zenostar® MO
Zirkonzahn Worldwide	-	Prettau® Anterior® Prettau® Zirkon ICE Zirkon Translucent
Zfx Zirkon	Zfx™ Zirconium	Zfx™ Zirconium allround Zfx™ Zirconium eff ect Zfx™ Zirconium Multilayer Zfx™ Zirconium BionX2

# Phosphate Cement for QM System

Single dose system acc. to MDR




An incorrect dosage demonstrably influences the mechanical and biological properties of all dental cements. Carboxylate cements react more strongly to too little liquid. Phosphate cements do not tolerate too little powder content. Even with a dosage spoon, deviations of 30% are not uncommon.

Single dosages guarantee standardization so that the same mixing ratio can always be achieved. In combination with the dexterity of the human hand, one and the same product can be dosed differently in order to enable flexible use – for the cementation of restorations, for core build-ups and for underfillings. READY2MIX fulfills documentation requirements for implant ID cards.


Advantages of **HOFFMANN'S** READY2MIX compared to manual mixing

- Always successful
- 33% shorter mixing time
- Optimum cement properties due to exact dosing
- Hygienic, moisture-proof powder packaging



Advantages of **HOFFMANN'S** READY2MIX over capsule or syringe systems

- One product for luting, lining and temporary filling
- No equipment required
- Economical over capsule and syringes
- 90% less waste
- Temperature and climate resistant





0%

BPA, monomeres,  
phenols, nanofillers

Compliant with quality management systems acc. to MDR

## READY2MIX – ZINC PHOSPHATE CEMENT NORMAL



## AREAS OF APPLICATION

- Lining for all filling materials (amalgam, composites, etc.)
- Cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal and ceramic
- Cementation of zirconium oxide, aluminium oxide and lithium disilicate ceramic
- Cementation of implant-supported crowns and bridges
- Cementation of orthodontic bands
- Cementation of retention pins and screws
- Core build-ups
- Long-term temporary fillings

## COMPOSITION

Powder: zinc oxide, magnesium oxide

Cement liquid: o-phosphoric acid

COMPARED TO  
CAPSULES & SYRINGES

- One system for lining, bonding and long-term provisional filling
- Individual consistency control
- No technical equipment needed
- 90% less waste

ADVANTAGES OF  
ZINC PHOSPHATE CEMENT

- No rubber dam
- Easy excess removal
- Low pulp toxicity
- Recommended for allergy sufferers
- Antibacteria effect

READY2MIX		
Order No.	Package size	Colour
80024-(colour)	Stick 20 x 1 g powder Tube 1 x 10 ml liquid	01, 03, 04

READY2MIX Pocket		
Order No.	Package size	Colour
8002355	Stick 5 x 1 g powder Tube 1 x 3 ml liquid	03



Zink polycarboxylate cement for luting and lining

**HOFFMANN'S CARBOXYLATE CEMENT**



**AREAS OF APPLICATION**

- Lining for all filling materials (amalgam, composites)
- Cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal and ceramic
- Cementation of orthodontic bands
- Cementation of retention pins and screws
- Core build-ups
- Long-term temporary fillings
- Fillings in deciduous teeth

**ADVANTAGES**

- No acid pain
- For sensitive tooth stumps
- Easy to use
- No etching, priming, bonding
- Compressive strength well above the material standard
- Chemical adhesion

**SETTING TIME FROM END OF MIXING:**



For cementation: 3 – 4 minutes  
For cavity linings: 2 – 3 minutes

**COMPOSITION**

Powder: zinc oxide, magnesium oxide  
Liquid: polyacrylic acid

CARBOXYLATE CEMENT		
Order No.	Package size	Colour
80312-(Colour)	1 x 100 g powder	01-15
8031300	1 x 40 ml liquid	---
80315-(Colour)	1 x 35 g powder 1 x 15 ml liquid	01-15





Compliant with quality management systems acc. to MDR

## AQUA CC - WATERMIX CARBOXYLATE CEMENT



### AREAS OF APPLICATION (as lining cement)

To compensate for polymerisation shrinkage under composite fillings. The special material properties of AquaCC make it a perfect lining material. AquaCC expands minimally and can thus compensate for shrinkage of composites and prevent the formation of marginal gaps. AquaCC also acts as a barrier against the free radicals formed during polymerization. AquaCC contains zinc and has an antibacterial effect.

### Further indications

- Cementation of restorations on implants
- Core build-ups
- Long-term temporary fillings

### ADVANTAGES

- No acid pain, especially suitable for sensitive tooth stumps and when working without anaesthetics
- Compressive strength well above the material standard
- Exact dosage and reproducible cement properties
- Ease handling
- Easy removal of excess material
- Chemical adhesion

### AESTHETICS / TIP

AquaCC is available in two colours: universal tooth-colour and light pink. The opacity of the lining cement ensures harmonious priming and retention of the colour stability even for longer periods of time. Due to the light pink hue, a particularly vivid effect can be achieved.

AQUA CC		
Order No.	Package size	Colour
8037010	Stick 10x1 g powder, pipette	universal
8037014	Stick 10x1 g powder, pipette	light pink

### COMPOSITION

Zinc oxide, magnesium oxide, polyacrylic acid



"Since bonding with Hoffmann's READY2PROTECT, No cases of periimplantitis occurred."

**DR. NICOLA MINOTTI**  
Dentist, Geneva

Compliant with quality management systems acc. to MDR

**CEMENT4IMPLANT - COPPERION CEMENT**



Bactericidal luting cement for implants

**AREAS OF APPLICATION**

- Permanent cementation of implant-supported crowns and bridges
- Cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal and ceramic
- Bonding of restorations made of zirconium oxide, aluminium oxide as well as lithium disilicate ceramic
- Cementation of orthodontic bands
- Cementation of retention pins and screws

**ADVANTAGES**

- Bacteria-proof closure
- Slow release of copper ions
- Pure micromechanical bonding
- If necessary, the restoration can be removed with an ultrasonic scaler
- X-ray visible for easy follow-up
- Easy excess removal
- Particularly well tolerated

**READY2PROTECT CEMENT4IMPLANT**

Order No.	Package size	Colour
8007010	Stick 10 x 1 g powder Tube 1 x 5 ml liquid	---

**COMPOSITION**

Powder: zinc oxide, magnesium oxide, copper salts  
Cement liquid: o-phosphoric acid



Phosphate cement with copper additive and bactericidal properties

## HOFFMANN'S COPPER CEMENT



### AREAS OF APPLICATION

- Fillings in deciduous teeth
- Lining for all filling materials (amalgam, composites)
- Cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal as well as metal ceramic and all-ceramic (zirconium oxide, aluminium oxide and lithium disilicate ceramic)
- Cementation of implant-supported crowns and bridges
- Cementation of orthodontic bands
- Cementation of retention pins and screws
- Core build-ups
- Long-term temporary fillings

### ADVANTAGES

- Dimensionally stable – no shrinkage and therefore no formation of marginal gaps
- Very easy removal of excess material
- Low risk of allergic reactions
- Bactericidal effect

COPPER CEMENT		
Order No.	Package size	Colour
8072203	1 x 100 g powder	03
8072300	1 x 40 ml liquid	---
8072403	1 x 30 g powder 1 x 13 ml liquid	03

### COMPOSITION

Powder: zinc oxide, magnesium oxide, copper salts  
Cement liquid: o-phosphoric acid

# Natural Healing Cement with Copperions

## Minimal invasive protocol

The protocol described below was developed by Dr. Jean-Pierre Eudier, Luxembourg, and Dr. Nicola Minotti, Switzerland. It describes a minimally invasive method for CP excavation and indirect pulp capping using zinc oxyphosphate cement with copperions.



1. Caries Profunda with initial (reversible) pulpitis (Fig. 1)

2. Atraumatic Preparation of the cavity (Fig. 2)

The cavity is carefully prepared, carious tissue is removed from the cavity walls. A layer of demineralized dentin is left on the pulp horn to avoid a pulp opening.

3. Application of a copper and copal barrier (Fig. 3)

Copal varnish mixed with a pinch of copperion cement powder is used as a liner. The ionized copal varnish is applied on the demineralised dentin. Copperions are transported into the dentin tubuli and the tubuli are sealed.

4. Mechanical Copperion Cement barrier (Fig. 4)

Lost dentin is replaced with a powder rich mixture of copperion cement. The cement can be used as a long term temporary filling. However a sustainable solution requires an enamel replacement with either gold or full ceramic. If budget restrictions prevent a sustainable restoration there are two other options: an amalgam rich in copper or a composite based filling (Preferably a classic etch & bond system without TEGDMA and HEMA).



Protocol video:  
[www.hoffmann-dental.com](http://www.hoffmann-dental.com)



„These cements are easy to use, economical and have a strong therapeutic potential. Nowadays they are unrivalled and very reliable.“

**DR. DANIEL BANDON**

(Lecturer, University of Marseille,  
Faculty of Paediatric Dentistry)

**0%**

BPA, monomers,  
phenols, nanofillers

## NATURAL HEALING CEMENT WITH COPAL AND COPPER

Compliant with quality management systems acc. to MDR

### READY2PROTECT COPPERION CEMENT + COPAL VARNISH



#### AREAS OF APPLICATION

- Healing cement to preserve vital teeth also with deep caries in combination with Hoffmann's Copal varnish using a copper barrier
- Minimal invasive caries treatment with modified ART technique
- Fillings in deciduous teeth
- Long-term temporary fillings, saliva-proof, with bactericidal effect
- Lining with long-lasting bactericidal effect for all types of teeth (also for vital teeth)
- Lining for all filling materials (composites and amalgam)
- Core build-ups

#### HOW IT WORKS

- The copal varnish infiltrates the dentin tubules and thus transports antibacterial copper ions into the demineralized dentin. The tooth pulp is supported in its primary function, so self-healing powers are activated.
- Missing dentin is replaced with READY2PROTECT. The cement layer forms a permanent barrier against secondary caries.

#### COMPOSITION

Powder: zinc oxide, magnesium oxide, copper salts

Cement liquid: o-phosphoric acid

Copal varnish: copal, ethanol

READY2PROTECT Pocket + Copal		
Order No.	Package size	Colour
8007020	Stick 11 x 1 g powder Tube 1 x 5 ml liquid Bottle 1 x 5 ml copal varnish	---





Pulp capping material with eugenol and propolis

## HOFFMANN'S PULPINE



### AREAS OF APPLICATION

- Caries profunda, indirect pulp capping
- Direct pulp capping, always without pulpitis

### HOW IT WORKS

The matrix with eugenol and propolis sets quickly and limits the possible side effects of alkaline calcium hydroxide, such as permanent stimulus or uncontrolled necrosis. The antibacterial properties of propolis (bee resin) supplement the disinfecting properties possessed by calcium hydroxide. Very good biocompatibility. Relieves pain.

### ADVANTAGES

- Biocompatible pulp capping alternative with propolis
- Controlled neutralization of calcium hydroxide
- Easy application
- Quick and complete setting
- Pain-relieving

### CONTRAINDICATION

Do not use if you are allergic to eugenol and/or bee products.

### SETTING TIME:

1:00 - 3:45 minutes

Hoffmann's PULPINE	
Order No.	Package size
84016	10 g powder
84017	10 ml liquid

### COMPOSITION

Powder: calcium compounds, zinc compounds

Liquid: ethanol, eugenol, propolis



Hydroxyapatite and propolis-containing pulp capping material and alternative to MTA

## HOFFMANN'S PULPINE MINERAL



### AREAS OF APPLICATION

- Caries profunda, indirect pulp capping
- Indirect pulp capping with reversible pulpitis
- Direct pulp capping, always without pulpitis
- Direct and indirect pulp capping after crown preparations

### HOW IT WORKS

Hydroxyapatite, which forms about 70% of natural dentin, is embedded in the composite of calcium hydroxide and propolis. It forms a very stable setting matrix without gaps and without necrosis formation. The antibacterial properties of propolis lead to a complication-free healing of infected pulp tissue.

### ADVANTAGES

- Antibacterial effect
- No necrosis formation
- Easy application
- Very good adhesion
- Quick setting

### CONTRAINDICATION

Do not use if you are allergic to eugenol and/or bee products.

### SETTING TIME:

1:00 - 3:30 minutes

Hoffmann's PULPINE MINERAL	
Order No.	Package size
84078	5 g powder
84077	10 ml liquid

### COMPOSITION

Powder: calcium compounds, hydroxyapatite

Liquid: ethanol, propolis



The alternative to MTA, biological bone cement on pure calcium phosphate and hydroxyapatite basis

## HOFFMANN'S ENDO REPAIR



### AREAS OF APPLICATION

- Vital extraction of deciduous teeth
- Perforation of the root canal
- Incomplete root growth, root filling
- Over-instrumentation of apex


### • ADVANTAGES

- Pure calcium phosphate matrix
- No necrosis formation
- Easy application
- Quick and complete setting

### HOW IT WORKS

The bone- or tooth-like matrix ensures maximum biocompatibility. Necrosis can be prevented and there is an almost asymptomatic healing of the pulp under the stable artificial bone layer.

**SETTING TIME:** 0:30 - 2:30 minutes



Hoffmann's ENDO REPAIR	
Order No.	Package size
84086	3,5 g powder
84087	10 ml liquid

### COMPOSITION

Powder: calcium phosphates, hydroxyapatite  
Liquid: distilled water (without preservatives)

# INDICATION TABLE

AREAS OF APPLICATION	PULPINE	PULPINE NE	PULPINE MINERAL	ENDO REPAIR	ENDO ABSORB	ENDO ABSORB PLUS	PROXI APEX	READY2PROTECT + Copal
Caries profunda [indirect pulp capping]	+++	+++	+++	++				++++
Indirect pulp capping with reversible pulpitis	++++	+++	+++					++++
Direct pulp capping [always without pulpitis]	+	++++	+++	+++				
Indirect pulp capping [after crown preparation]	++	+++	++++	+++				
Direct pulp capping [after crown preparation]		+++	++++	+++				
Exsudate formation during root canal treatment					+++	++++		
Bleeding during root canal treatment					+++	+++		
Root filling [vital extirpation]							++	
Final root filling [after gangrene treatment]							++	
Vital extraction of deciduous teeth				++++				
Perforation of the root canal [temporary]				++++				
Perforation of the root canal [final]				++++				
Incomplete root growth [root filling, temporary]				++++				
Incomplete root growth [root filling, final]				++++				
Over-instrumented apex				++++				

+ possible  
 ++ good  
 +++ very good  
 ++++ excellent

## Bentonit - „Healing Clay“



Temporary root canal therapeutic for absorbing the exudate. Radiopaque

### HOFFMANN'S ENDO ABSORB PLUS



#### AREAS OF APPLICATION


- Treatment of infected root canals after pulp removal or gangrene treatment
- Treatment of exudate formation and bleeding after root canal treatment

#### HOW IT WORKS

Temporary root canal therapeutic with antimicrobial effect and a unique feature: It absorbs exudate up to an amount of 120% of its own weight and over a period of 3 weeks. Contains calcium hydroxide. Radiopaque.

#### ADVANTAGES

- Absorption of the exudate up to 120% of its own weight
- Antimicrobial
- Very good compatibility
- Absorption capacity of 3 weeks
- Radiopaque



**SETTING TIME:** 1:00 - 4:00 minutes

Hoffmann's ENDO ABSORB PLUS	
Order No.	Package size
84040	5 g powder
84039	10 ml liquid

#### COMPOSITION

Powder: calcium compounds, bentonite, radiocontrast agent  
Liquid: alcool





Permanent root canal filling material based on zinc oxide and eugenol with antibacterial effect. Radiopaque

## HOFFMANN'S PROXI APEX



### AREAS OF APPLICATION

Permanent root canal filling material

- After vital extirpation
- After gangrene treatment
- Before apicoectomy

### HOW IT WORKS

Hoffmann's PROXIAPEX is a eugenol-containing root filling material for permanent treatments. During the 20 minute setting process the pH value is neutralized. The material exhibits a mild antibacterial effect and is radiopaque.

### ADVANTAGES

- Very good tolerability
- Antibiotic and cortisone-free mild antibacterial effect
- Neutral pH value
- No irritation of the apical soft tissue after setting
- Good canal wall resistance
- Tight sealing of the dentin canals thanks to increased viscosity
- Quick setting
- High radiopacity

### CONTRAINDICATION

Do not use if you are allergic to eugenol.

SETTING TIME: 4:00 - 15:00 minutes

Hoffmann's PROXI APEX	
Order No.	Package size
84026	10 g powder
84027	10 ml liquid

### COMPOSITION

Powder: zinc compounds, calcium compounds, zirconium oxide  
Liquid: eugenol



Self-cure sealing compound for the temporary restoration of cavities.



**HOFFMANN'S EASYFILL AKUT**



**AREAS OF APPLICATION**

- Ready-to-use temporary seal (hydrophilic) for cavities
- Emergency tooth filling in the context of the dental emergency service
- secure seal for inlay medications

**HOW IT WORKS**

EASY FILL akut is a ready-to-use temporary filling paste with a neutral pH value belonging to the class of zinc sulphate cements. The temporary sealing compound is taken out of the jar with an instrument and placed directly into the cavity. EASY FILL akut can be used under specific situations for instance as an emergency tooth filling in the context of the dental emergency service. Thus preventing food and bacteria from getting inside the tooth cavity and leading to further oral problems. EASY FILL akut does not need to be heated or made ready for use by other means. EASY FILL acute is a paste that cures by itself when moisture (saliva) is present.

**ADVANTAGES**

- secure marginal seal
- eugenol-free
- PH neutral
- odourless
- self-curing
- radiopaque
- easy to remove

Hoffmann's EasyFill akut	
Order No.	Package size
8100400	40 g paste

**WARNINGS / CONTRAINDICATIONS**

- Do not apply directly to the opened pulp.
- Cover dentin close to the pulp beforehand, e.g. with Hoffmann's Pulpine/PULPINE NE
- Should only be used after prior risk/benefit assessment in patients with severe renal dysfunction (nephritis and/or nephrosis).

**FURTHER INFORMATION FOR PATIENT SAFETY EDUCATION**

Water can be sipped in small portions after treatment. Curing time is 30 minutes, nonetheless chewing should be avoided for 2 hours after treatment.

The exact life span of this temporary filling may differ from one person to the other. However, EASY FILL akut should not stay in the cavity for longer than 30 days.

Patients wit severe renal dysfunction (nephritis and/or nephrosis) should only be treated after prior detailed consultation.

**COMPOSITION**

Zinc oxide, Potassium sulphate, Zinc sulphate, Vinyl chloride vinyl acetate



Traditional ZOE Cement

## HOFFMANN'S ZINKOXID EUGENOL SET



### AREAS OF APPLICATION

- Temporary Obturation
- Temporary Root Canal Seal
- Intermediate restorative material

### HOW IT WORKS

Zinc oxide and eugenol (natural clove oil) for the preparation of a classic zinc oxide eugenol paste (cement). (ZOE) Both components are mixed to form a Zinc-Oxide-Eugenol Complex/ Cement.

Setting and mixing time depend on powder to liquid ratio. The setting time (tested according to ISO 3107) is 2h30 - 8:00h according to the recommended Powder/Liquid ratio.

### FURTHER INFORMATION FOR PATIENT SAFETY EDUCATION

The strength of the setting cement increases slowly. Water can be sipped in small portions after treatment. Chewing should be avoided depending on the powder/liquid ratio selected, i.e. 2.5 to 8 hours after treatment. Hoffmann's Zinoxid Eugenol should not be used if the patient is known to be allergic to eugenol. In case of doubt, an epicutaneous or LTT test should be carried out before the planned treatment.

### ADVANTAGES

The traditional formula (without additives) is characterised by an exceptionally long curing time, delivering an initial soothing effect thanks to the innate analgesic and antibacterial properties of the two components.

### WARNINGS / CONTRAINDICATIONS

Do not apply directly to the opened pulp.  
Cover dentin close to the pulp beforehand, e.g. with Hoffmann's Pulpine/Pulpine NE

Eugenol impairs the setting of (artificial) resin cements/composites and can also soften composites that have already set.

- Before composite fillings can be applied,
- 1) a lining material must be placed (e.g. sandwich-technique with Hoffmann's Cement) and
  - 2) residual eugenol must be neutralised. For example, by means of the total-etch technique or by treatment with ozone.

Do not use after the expiry date.

Hoffmann's Zinkoxid Eugenol Set	
Order No.	Package size
8001002	30 g powder, 15 ml liquid

### COMPOSITION

Powder: Zinc Oxide (Ph. Eur.)

Liquid: Eugenol (natural clove oil) Ph. Eur.



Temporary sealing material with antibacterial properties

HOFFMANN'S GUTTAPERCHA



AREAS OF APPLICATION

- Temporary filling of cavities
- Sealing of holes in implant screws during the healing phase

ADVANTAGES

- Genuine raw gutta-percha
- High biocompatibility
- Antibacterial effect
- Easy to remove without residues
- Easy removal from the cavity with a probe or an excavator
- Eugenol-free

GUTTAPERCHA		
Order No.	Package size	Colour
82303	Sticks, 110 g	---

COMPOSITION

Zinc oxide, gutta-percha, paraffin, beeswax



Temporary filling cement with antibacterial properties

## HOFFMANN'S PROVISIONAL CEMENT



### AREAS OF APPLICATION

- Classic short-term temporary closure of cavities

### ADVANTAGES

- Easily removable from the cavity
- Saliva-proof over a period of one month
- Biocompatible due to exclusively inorganic ingredients

For classic short-term temporary fillings, Hoffmann's Provisional Cement is saliva-proof over a period of one month and can be used for temporary sealing. The durability of the temporary seal strongly depends on the oral flora and oral hygiene, as well as the diet of the individual patient.

**SETTING TIME:** 1:00 - 3:00 minutes

TEMPORARY CEMENT		
Order No.	Package size	Colour
8100203	1 x 50 g powder	---
8100300	1 x 40 ml liquid	---

### COMPOSITION

Powder: zinc oxide, aluminium hydroxide

Liquid: o-phosphoric acid





Stents thermoplastic impression compound for precise functional impressions

**HOFFMANN'S IMPRESSION COMPOUND RED - STENTS**



**AREAS OF APPLICATION**

- Individual impressions
- Bite registration, also of edentulous mandibles before the production of special splints

Also suitable for small and highly sophisticated modelling techniques and modelling techniques used in jewellery workshops.

**PROCESSING TEMPERATURE**

58 – 62 °C

**ADVANTAGES**

- Easy to heat
- Easy to knead
- Does not stick to gloves, teeth or stumps
- Bite registration can be easily repositioned, corrected and milled
- Dimension-true reproduction of detail
- Odourless
- Composition of natural raw materials with a high degree of pharmaceutical-grade purity
- Physiologically harmless
- Compostable

IMPRESSION COMPOUND RED		
Order No.	Package size	Colour
82604	6 plates	red

**COMPOSITION**

Stearin, copal, talc



Stents thermoplastic impression compound for precise impressions and occlusal bite registration

HOFFMANN'S IMPRESSION COMPOUND GREEN - STENTS



AREAS OF APPLICATION

- Correction of impressions
- Enhancing the impression edges
- Occlusal bite registration
- Copper ring impressions

PROCESSING TEMPERATURE

58 – 62 °C

ADVANTAGES

- Easy to heat
- Excellent flowability
- Even deformability
- Bite registration can be easily repositioned, corrected and milled
- Rapidly hardens in cold water
- Composition of natural raw materials with a high degree of pharmaceutical-grade purity
- Physiologically harmless
- Compostable

IMPRESSION COMPOUND GREEN		
Order No.	Package size	Colour
82607	15 sticks	green

COMPOSITION

Stearin, copal, talc, carnauba wax



## HOFFMANN'S CASTING WAX - INLAY WAX (BLUE OR IVORY)



### AREAS OF APPLICATION

- Stick casting wax with excellent natural plasticity, excellently suited for the modelling of inlays, onlays, crowns and bridges
- In many cases also suited for small and complicated modelling techniques used in jewellery workshops

### MELTING POINT

ca. 75 °C

BLUE - traditionally for cast gold fillings

IVORY - for ceramic and synthetic restorations

### ADVANTAGES

- Particularly suitable for the traditional direct wax technique
- Dimensionally stable
- Residue-free combustion
- Composition of natural raw materials with a high degree of pharmaceutical-grade purity
- Physiologically harmless
- Biodegradable

### TIP

Processing is carried out by means of a Bunsen burner, open flame and electrical wax knife.

CASTING WAX		
Order No.	Package size	Colour
83001	20 sticks, 30 g	blue
83101	20 sticks, 30 g	ivory

### COMPOSITION

Paraffin, carnauba wax, dammar resin, stearin



## HOFFMANN'S COPAL VARNISH - CAVITY VARNISH



### AREAS OF APPLICATION

- Disinfecting cavity protection varnish for the avoidance of post-operative sensitivity
- Cavity protection varnish before cement fillings
- Sealing of dentinal tubules
- Lining of cavities before amalgam fillings and restorations made of gold to protect against corrosion products
- Insulation against thermal irritation
- Protective varnish for the surfaces of cement fillings and plaster models

### ADVANTAGES

- Made of natural raw materials
- Ready-to-use solution
- Contains no chloroform or ethyl ether
- Fluoride free

COPAL VARNISH		
Order No.	Package size	Colour
81901	50 ml, liquid	---

### COMPOSITION

Copal, ethanol



Adhesive wax speciality – breaks like glass

**HOFFMANN'S ADHESIVE WAX (RED OR YELLOW)**



**AREAS OF APPLICATION**

- Blocking out undercuts on the model
- Fixation of parts to be soldered before making the soldering model
- Crown and bridge modellation: Fixation of model to a casting funnel, blocking and securing the cap to a bridge
- Re-fixation of severed bridge modellations in order to prevent tension in the wax
- Fracture repairs of prostheses: Fixation of broken denture bases or chipped teeth for subsequent making of model or matrix
- Fixation of the maxilla and mandible for insertion in an articulator

**MELTING POINT**

ca. 70 °C

**ADVANTAGES**

- Particularly suitably for the traditional direct wax technique
- Stick shape very easily workable by heating the wax directly over a flame (drips) or controlled melting on a hot instrument
- Excellent adhesive properties; adheres to gypsum or other dental materials
- Made of pure natural raw materials
- 100% error control, no undetected warping of modellations or fixations: the material is hard and brittle at room temperature; the wax does not bend but breaks under load
- Easy removal using steam
- Residue-free combustion – with the pleasant smell of beeswax

ADHESIVE WAX		
Order No.	Package size	Colour
82914	14 sticks, 70 g	red
82913	14 sticks, 70 g	yellow

**COMPOSITION**

Colophony, bees wax, carnauba wax



Allergy-free and sustainable work in a dental laboratory with natural shellac base plates

## HOFFMANN'S SHELLAC BASE PLATES - PINK



### AREAS OF APPLICATION

- Individual impression trays
- Making dentures: Base material for bite impressions and tooth lineups

### FEATURES

- Easily and uniformly processed with heat
- Dimensionally stable (retains shape during try-in in patient mouth)
- Available in two versions for maxilla and mandible
- Material thickness approx. 1.3 mm

### Advantages compared to light-curing and autopolymerizing materials

- Corrective measures can be taken without material-loss through re-heating
- Natural raw materials with a high degree of pharmaceutical-grade purity
- Physiologically harmless
- Biodegradable
- Skin-friendly – no known allergies
- Pleasant odour

### PROCESSING RECOMMENDATIONS

Hoffmann's shellac base plates should be heated very gently with a customary CE certified hot air blower. Hot air provides a uniform spreading of heat and thus reduces the risk of burning.

### COMPOSITION


Shellac, stearin, talc, mica



Video on the subject:  
[www.hoffmann-dental.com](http://www.hoffmann-dental.com)

SHELLAC BASE PLATES	
Order No.	Package size
82824	8 plates maxilla + 4 plates mandible
82825	12 plates maxilla
82826	12 plates mandible
8282104	50 plates maxilla + 50 plates mandible
8282105	100 plates maxilla
8282106	100 plates mandible





Recommended by:  
**agkeramik**

In clinical use since 35 years. (Roulet 1989 Theiss 2009)

**HOFFMANN'S SILANE - COUPLING AGENT**



**AREAS OF APPLICATION**

- To increase the chemical bonding force during adhesive cementation of silicate ceramics, silicated oxide ceramic or metal restorations with methacrylate bearing cementation materials (composites, compomers and ormocers)
- As a ceramic-composite bonding agent for ceramic repairs

**Benefits of Hoffmann's two-component silane**

- Full reaction force as the reagents are mixed freshly before each use
- Longer shelf life after opening the bottle
- Optimal chemical bonding surface with simultaneous good surface wettability
- Prevention of adhesion defects

SILAN		
Order No.	Package size	Colour
82101	Component A, 5 ml Component B, 5 ml	---

**COMPOSITION**

Acetic acid in ethanol solution, 3-methacryloyloxypropyltrimethoxysilane in ethanol solution

On hot summer days  
put a cool pack under the glass plate!

TIPP



HOFFMANN'S EXPERT MIXING TOOLS

TIPP

Glass mixing block – solid crystal block

Shatter-proof, resistant to abrasion even from metal spatulas, and easy to clean. The perfect basics for mixing all Hoffmann's cements.

Throughout the entire mixing process the glass maintains a constant temperature even in the case of heat from reaction. The glass block provides plenty of space for spreading and allows safe and ergonomic handling.

Cement spatula – heavy version

The cement spatula provides high work comfort and facilitates the mixing of cements – even with very solid, packable lining consistencies.

EXPERT MIXING TOOLS		
Order No.	Tool	Dimensions and weight
82900	Mixing block made of crystal glass, solid	150 x 75 x 20 mm, 596 g
82908	Cement spatula, heavy version	approx. 172 mm, approx. 23 g

## HOFFMANN'S Backlist



**Hoffmann's UNIVERSAL CEMENT**



100 g Powder  
REF 80032-[colour]

Farbnr. 01-15



**Hoffmann's PULPINE NE**



10 g Powder  
REF 84056



10 ml Liquid  
REF 84037



**Hoffmann's ENDO ABSORB**



5 g Powder  
REF 84035



10 ml Liquid  
REF 84057

## HOFFMANN'S Gadgets



Hoffmann's  
DENTAL MEMO  
Memory game



Hoffmann's  
CROWN & BEARD



Tales of Hoffmann's  
„The Inventor“  
DE, FR, IT, RU, ES, NL



Hoffmann's SET  
OF POST CARDS  
(7 motifs) with display



Tales of Hoffmann's  
„In Africa“  
DE, FR

# Advisors / Accountants

DR. JEAN-PIERRE EUDIER

*Dentist & Senior Advisor*  
FR, EN

54 Rue de la Montagne  
L-3259 Bettembourg, Luxembourg  
T +352 691 18 34 12  
eudier@hoffmann-dental.com

ANTONY FILHOL

*International Key Account*  
DE, EN, ES

Komturstr. 58-62  
D-12099 Berlin, Germany  
T +49 30 - 82 00 99 - 0  
a.filhol@hoffmann-dental.com

## HOFFMANN'S Nomenclature



Powder stick pack



Copper ion



Eugenol



Sticks



Liquid tube



Radiopaque



Propolis



Round plates



Powder glass



Set



Syringe



Lower jaw



Liquid bottle



Luting



Ozone therapy



Upper jaw



Dosing bottle



Lining



Algae



Component



Paste-like



Colours of cement

# Order Form



**HOFFMANN'S**

☐ New customer ☐ Customer ID \_\_\_\_\_

☐ Dentist ☐ Dental technical ☐ Dental assistant ☐ College ☐ Dealer

Clinic / Company name:			
Contact person:	<input type="checkbox"/> Mrs	<input type="checkbox"/> Mr	
Street / N°:			
Postcode / Country / City:			
Tel.:		<input type="checkbox"/> EU	<input type="checkbox"/> other country
Fax:			
Email:			
Depot:			

**www.hoffmann-dental.shop**

VAT No.: \_\_\_\_\_

Order No.	PRODUCT	COLOUR	QUANTITY	PRICE IN EURO [€]

☐ Yes, I would like to receive the news letter by e-mail.

Total

Order fax: +49 30 82 28 14 5

E-Mail:  
sales@hoffmann-dental.com

Hereby I confirm my order.

Date \_\_\_\_\_ Signature \_\_\_\_\_

## Notes



# Contact

For detailed information on materials, processing and sources of supply, please contact our medical device consultants & our sales team.

## POSTAL ADDRESS:

Hoffmann Dental Manufaktur GmbH  
Wangenheimstraße 37/39  
14193 Berlin  
Germany  
Tel.: +49 30 82 00 99 0  
Fax: +49 30 82 28 145  
E-Mail: [sales@hoffmann-dental.com](mailto:sales@hoffmann-dental.com)

**[www.hoffmann-dental.com](http://www.hoffmann-dental.com)**

## IMPRINT

© 2023 Hoffmann Dental Manufaktur GmbH  
Wangenheimstraße 37/39  
14193 Berlin | Germany  
[www.hoffmann-dental.com](http://www.hoffmann-dental.com)

Publisher: Yvonne Hoffmann

Illustration: Norbert Folberger (pp. 9, 10, 12, 21)

Reprint, even in part, only allowed with permission of Hoffmann Dental Manufaktur GmbH.

Printed in Germany



**FACTORY SALES:**  
Konturstraße 58-62  
12099 Berlin

By appointment only.

 Otto Hoffmann

 HoffmannDental

 Hoffmann Dental

 hoffmann\_dental

#hoffmannndental