



# HOFFMANN'S UNIVERSAL CEMENT



Zink Phosphat Cement

---

## Gebrauchsanweisung für den Zahnarzt

Zusammensetzung: Pulver enthält: Zinkoxid, Magnesiumoxid

Nur zur Verwendung mit den Flüssigkeiten:

Hoffmann's Cement normalhärtend

Hoffmann's Cement schnellhärtend

Hoffmann's Carboxylat Cement

Anwendungsgebiete:

Befestigung von Inlays, Kronen und Brücken (auch aus Zirkonoxid);

Befestigung von orthodontischen Bändern;

Einzementierung von Wurzelstiften und Schrauben;

Unterfüllungen und temporäre Füllungen (nicht mit der Flüssigkeit Hoffmann's Cement schnellhärtend)

**BITTE BEACHTEN SIE DIE GEBRAUCHSANWEISUNG DER VERWENDETEN FLÜSSIGKEIT!**



# HOFFMANN'S UNIVERSAL CEMENT



Zinc Phosphate Cement

---

## Instruction for Use for the dentist

Compositions: powder contains: zinc oxide, magnesium oxide

Only to use with the liquids:

Hoffmann's Cement normal setting

Hoffmann's Cement quick setting

Hoffmann's Carboxylate Cement

Areas of application:

Final cementation of cast alloy inlays, crowns and bridges; cementation of zirconia ceramic crowns and bridges; cementation of orthodontic bands, retention pins and screws; cavity lining as base for all filling materials; temporary fillings even for longer terms (not applicable with Hoffmann's Cement quick setting liquid).

**PLEASE NOTE THE INSTRUCTION FOR USE OF THE USED LIQUID!**